

The GaBP Ring makes the difference.

Banded gastric bypass with Fobi technique



The GaBP Ring makes the difference.


Banded gastric bypass with Fobi technique




GaBP Ring

5.5 - 7.0cm circumference

- pseudopyloris of proximal pouch
- controls emptying of the pouch
- enables long-term effectiveness of gastric restriction






Banded gastric bypass with Fobi technique


The major differences:

- standardized shape (4 sizes):
reproducible results
- non-elastic material:
effective restriction without erosion
- autolock system:
reliable closing of GaBP Ring


1 guiding introducer through lesser omentum



2 align both ends of GaBP Ring for insertion



3 sutured GaBP Ring in place



source: Fobi article



- 4 GaBP Ring Autolock-Kit
Gastrostomy Site Marker
- 5 Lock Mechanism
- 6 GaBP Ring in closed position




The GaBP Ring makes the difference.

Banded gastric bypass with Fobi technique



Long-Term Outcomes after Gastric Bypass

Graph showing median weight of 342 patients before and after gastric bypass.



source: Stubbs article Presented 9th World Congress IFSO 8-11; 2004 · Obes Surg 2005; 15: 155-163

distributed by



Bariatric Solutions GmbH
Kaltenbacherstrasse 24
CH - 8260 Stein am Rhein

Phone: +41 - 52 741 - 65 00
Fax: +41 - 52 741 - 65 02
info@bariatric-solutions.com
www.bariatric-solutions.com